

# Comparative study of Exhaust Air Dust (EAD®) testing and traditional bedding sentinels for health monitoring in IVC rack systems

Dr. S. Durand\*, Dr. M. Brielmeier \*\*, Dr. A. Gobbi\*\*\*, Dr. K. Henderson\*\*\*\*

\*Veterinary and Professional Services, Charles River, Lyon, France

\*\*Research Unit Comparative Medicine, Neuherberg, Germany

\*\*\*Cogentech, Milano, Italy

\*\*\*\*RADS, Charles River, Wilmington, USA

charles river

## 1 INTRODUCTION

Alternative health monitoring such as Exhaust Air Dust (EAD®) PCR testing is a new method designed to improve the detection of pathogens. EAD® can be easily collected with in-line capture media available in different IVC rack systems.

Three independent research groups (Max Planck Institute-Frankfurt; Helmholtz Zentrum- Munich, Cogentech Research Center- Milano) performed a health monitoring comparison study.

They compared the accuracy of environmental health monitoring by using an in line EAD® filter in IVC racks (N=26) to the traditional bedding sentinel methodology.

## 2 MATERIAL AND METHODS

Each facility cleaned between 2-5 air handling units (AHU, Tecniplast) with different types of rack clusters (between 1 to 4 single or double racks linked to the AHU). In total, the study included 9 AHU's with 26 IVC racks. Each sample type was pre-processed and DNA/RNA was extracted using a total nucleic acid isolation kit. Sample nucleic acid was screened by individual single-plex PCR assays on a real time open array PCR platform. Semi quantification was performed to define copy number per reaction. Test specimens at 1, 2 and 3 months included base line swabs, prefilter swabs (N=7) and in-line EAD® filters (Interceptor, N=15) for environmental samples. At 3 months, soiled bedding sentinels (N=26) were sent to laboratories for traditional testing including bacteriology, serology, parasitology and PCR (Fig. 1).

## 3 RESULTS

Post cleaning, baseline PCR tests were all negative except for one sample which was found positive for *Pasteurella pneumotropica* at very low copy number (1-10 copies)(data not shown).

Table 1: Percentage of detected agents

Agents	Soiled bedding sentinels	Interceptor (N=15)	Prefilter swabs (N=7)
Virus			
MNV	42,3%	46,7%	42,86%
bacteria			
<i>Pasteurella pneumotropica</i> -Jawetz	3,8%	53,3%	85,71%
<i>Pasteurella pneumotropica</i> -Heyl	0,0%	60,0%	57,1%
<i>Helicobacter</i> genus	23,1%	100,0%	100,0%
<i>Helicobacter hepaticus</i>	7,7%	60,0%	85,7%
<i>Helicobacter ganmani</i>	0,0%	80,0%	100,0%
<i>Helicobacter mastomyrinus</i>	0,0%	46,7%	57,1%
<i>Helicobacter typhlonius</i>	0,0%	73,3%	100,0%
<i>Pseudomonas aeruginosa</i>	7,7%	0,0%	42,9%
<i>Staphylococcus aureus</i>	0,0%	20,0%	28,6%
parasites			
<i>Myocoptes musculus</i>	0,0%	6,7%	14,29%
<i>Entamoeba</i>	3,8%	40,0%	28,57%
Average	7,4%	48,9%	61,9%

Bad Good

The data showed a significant increase in detection rate with the filter media. Overall, the average percentage of detected agents using the media filter was 7 time higher as compared to the sentinels (Table 1).

Most of the agents were detected after 1 month of exposure (Fig. 2). Agents such as *S. aureus*, *P. aeruginosa* showed a delayed detection at 2 or 3 months (Fig.2).

*Helicobacter* spp showed the highest copy number where as *S. aureus* showed the lowest copy number (Fig.3). Opportunistic agents such as *S. aureus* shedded at a very low rate in immunocompetent mice.

Figure 1: Sample collection




Figure 2: Percentage of detected agents





Figure 3: Average copy number/µl of nucleic acid



## 4 CONCLUSION

- EAD® testing is an effective method for health monitoring: it improves the detection of pathogens as compared to soiled bedding sentinels in IVC racks
- EAD® testing can help reduce the number of sentinel animals needed per experiment: It's a positive way to apply the 3R's in an animal facility.

### Acknowledgements :

A. Yarali, Dr. Mombaerts from Max Planck Institute, Frankfurt  
 Dr Miller from Helmholtz Zentrum, Munich  
 DMV G. Milite from Techniplast  
 C. Perkins, the molecular diagnostic team, M. Sender, M. Galasso from CRL