

in vitro PDX models: 3D cultured patient-derived tumors for compound evaluation

AACR 2017
5770

Sander Basten¹, Bram Herpers¹, Julia Schueler², Torsten Gieseemann², Leo Price¹
¹OcellO B.V., Leiden, The Netherlands. ²Charles River Laboratories, Freiburg, Germany


Sander.Basten@OcellO.nl

Background: Patient-derived xenograft (PDX) models in immune-compromised mice allow propagation of and compound testing in human-derived tumors *in vivo*. To expand the potential of these human-relevant PDX models, we sought to develop 3D *in vitro* culture methods for PDX-derived tumor cells that show *in vivo*-like growth characteristics, invasion and responses to therapeutics. In combination with advanced 3D image analysis methods, we created a unique high throughput *in vitro* PDX screening platform that not only allows efficient identification of active and selective molecules but also enables selection of the optimal PDX tumor models for subsequent validation of candidates *in vivo*.


Results: Each PDX model has its own unique growth characteristics. Hydrogel and growth media composition were optimized to support growth of tumor tissues *in vitro* from cells derived from bladder, stomach, breast, pancreas, colon and lung cancer PDX tumors. Tumor tissues were cultured in a 384-well format and used to test chemotherapeutics (e.g. 5-FU, doxorubicin, paclitaxel, cisplatin), small molecules (e.g. erlotinib, lapatinib, trametinib, everolimus), antibodies (e.g. cetuximab, trastuzumab) and antibody-drug-conjugate (ADC, T-DM1) dose ranges. Using OcellO's 3D image analysis platform, Ominer, tumoroid growth, cell proliferation, apoptosis, invasion, cell polarity, differentiation and other aspects of cell and tissue architecture were analyzed and the effects of compound exposure on tumoroids was determined. By performing feature training based on reference compounds, we selected ±10 morphological features (out of more than 500) to generate a phenotypic signature that described the unique phenotypic change induced by each compound. Different compounds that target the same molecule were found to induce a similar morphological change whereas compounds with off-target effects could be discriminated. This approach enabled a high resolution evaluation and comparison of compound activity in an automated manner.

Conclusions: We established several PDX model-derived 3D tumor cultures in which standard-of-care and novel therapeutic agents (small molecules, antibodies and ADCs) can efficiently be screened, based on therapeutically relevant parameters and their changing morphological profile. This method enables both the *in vitro* selection of promising compounds in a pre-clinically relevant setting and the selection of optimum PDX tumor models for follow-up *in vivo* studies. This highly translational *in vitro-in vivo* PDX pipeline is expected to reduce attrition and increase efficiency in early drug-discovery.

Ex vivo propagation of human-relevant tumoroids



3D Tumoroid growth in extracellular matrix rich hydrogels




In vitro growth phenotypes of tumoroids cultured from murine PDX of various tissue origin and histotypes. Each model is grown in a optimized mix of extracellular matrix and growth factors. Retention of a relevant tumor marker (HER2) in *in vitro* tumoroids is shown (IF).

Developed in vitro PDX

ID	Tumor Type	Model Number
1	Bladder (AC)	BYE 138
2	Bladder (AC)	BYE 1029
3	Stomach (Adeno (AC))	GXA 3068
4	Stomach (Adeno (AC))	GXA 3069
5	Stomach (Adeno (AC))	GXA 3054
6	Stomach (Adeno (AC))	GXA 3067
7	Stomach (AC)	GXF 281
8	Lung (NSCLC, AC)	LXFA 677
9	Lung (SCC)	LXFE 211
10	Lung (AC)	LXFE 1029
11	Breast (HER2-)	MAXF 1169
12	Breast (HER2-)	MAXF 1329
13	Breast (TNBC)	MAXF 1365
14	Breast (HER2+)	MAXF 1389
15	Breast (TNBC)	MAXF 1011
16	Colon (AC)	OXF 1259
17	Colon (AC)	OXF 1733
18	Colon (AC)	OXF 1738
19	Colon (AC)	OXF 2061
20	Colon (AC)	OXF 2066
21	Colon (AC)	OXF 2093
22	Colon (AC)	OXF 2058
23	Colon (AC)	OXF 1103
24	Colon (AC)	OXF 699
25	Colon (AC)	OXF 2093
26	Pancreas (AC)	PAXF 1870
27	Pancreas (AC)	PAXF 1869
28	Pancreas (AC)	PAXF 1900
29	Pancreas (AC)	PAXF 2005
30	Pancreas (AC)	PAXF 2198


Image analysis and segmentation



Nuclei and actin signals are segmented using OcellO's analysis platform, Ominer. This allows for multiparametric measurement of morphological features to describe the change in the complete 3D structure, including:


- Number of nuclei, shape, distribution
- Fraction apoptotic nuclei
- Tumoroid volume, roundness, lumen content, cell polarity
- Microtissue branching and invasion.

After segmentation: measurement of relevant features




Dose dependent effects for two selected features (>500 are measured) in tumoroids derived from mammary model MAXF-1384 (TNBC). Depicted are end-points of a 7 day exposure to a panel of chemotherapeutic compounds and small molecules.

Differential drug responses



Model preselection; calculated IC50 values for reference compounds are compared amongst PDX models, identifying compound-sensitive models.

Morphologic clustering of compounds



Compound classes group in discriminatory clusters depending on the mechanism of action, allowing characterization of novel compounds by morphology.

Conclusions

Our platform allows drug sensitivity screening of Patient-derived tumor material for the purpose of pre-clinical drug development. In collaboration with Charles River, OcellO offers companion *in vitro - in vivo* PDX, combining human relevant tumor models with high-content compound evaluation.

- High-content pipeline
- Well-characterized PDX models
- Multiparametric analysis for selected features
- Amendable to small molecules, antibodies, ADC, etc
- Translational follow-up in **matching** *in vivo* models