

RABBIT ANTISERUM

Rabbit p27 antiserum was prepared in rabbits using purified virion protein (m.w.=27,000) extracted from avian myeloblastosis virus. Antiserum reacts with all subgroups of Leukosis virus.

Rabbit p27 Antiserum

LiquidMaterial No. 10100766

LyophilizedMaterial No. 10100763I

These reagents are derived from our Rabbit p27 antiserum. Normal useful working dilution of Anti p27 IgG and Anti p27 conjugate are approximately 1:1000.

Rabbit Anti-p27 IgGMaterial No. 10100771 (Lyophilized)

Rabbit Anti-p27 IgG Conjugated to HRPMaterial No. 10100770 (Lyophilized)

p27 Positive TC Control.....Material No. 10100764 (Lyophilized)

p27 Negative TC ControlMaterial No. 10100765 (Lyophilized)

Charles River | Avian Vaccine Services

Franklin Commons, 106 Route 32, North Franklin, CT 06254, USA

Phone: 1-800-772-3271 (USA toll-free) • 1-860-889-1389 • Fax: 1-860-889-1991 • E-mail: csdspf@crl.com

Wire Transfer: JP Morgan, New York, NY, USA, DDA Account #799-761499, ABA# 021000021, SWIFT #CHASUS33

Checks: Please mail to Charles River Laboratories, G.P.O. Box 27812, New York, NY 10087-7812