

MAIN OFFICE  
 Franklin Commons  
 106 Route 32  
 North Franklin, CT 06254  
 Tel: (860) 889-1389  
 Fax: (860) 889-1991



LABORATORY  
 67 Baxter Road  
 Storrs, CT 06268  
 Tel: (860) 429-7261  
 Fax: (860) 429-2572

## SPF RESEARCH EGGS

### QUALITY CONTROL SHEET

Flock \_\_\_\_\_ Hatch Date \_\_\_\_\_ Country of Origin \_\_\_\_\_

This chicken flock has tested Negative to the following agents by the test methods indicated below.

Agent	Antigen	Test
Avian Adenovirus Group I	CELO-Phelps	AGP
Avian Adenovirus Group II (HEV)	Domermuth	MFIA
Avian Adenovirus Group III (EDS)	CLKK115D	HI
Avian Encephalomyelitis	van Roekel	MFIA
Avian Influenza (Type A)	T/W/66	AGP, ELISA
Avian Nephritis Virus	G4260	MFIA
Avian Paramyxovirus Type 2	Yucaipa	MFIA
Avian Reovirus	S 1133	MFIA, AGP, IFA
Avian Rhinotracheitis Virus	UK	ELISA
Avian Rotavirus	Ch-2	AGP
Avian Tuberculosis	M. avium	CO, PM
Chick Anemia Virus	DelRose	IFA
Endogenous GS Antigen	p27	ELISA
Fowl Pox	Conn	MFIA, CO
Hemophilus paragallinarum	Serovars A,B,C	CO
Infectious Bronchitis - Ark.	Ark 99	MFIA
Infectious Bronchitis - Conn.	Conn A5968	MFIA
Infectious Bronchitis - JMK	JMK	MFIA
Infectious Bronchitis - Mass.	Mass 66579	MFIA
Infectious Bursal Disease Type 1	M4040(2512)	MFIA, AGP
Infectious Bursal Disease Type 2	M4040(2512)	AGP
Infectious Laryngotracheitis	UC A92430	AGP, ELISA
Lymphoid Leukosis A, B	RSV-RAV A, B	MFIA
Avian Lymphoid Leukosis Virus J (ALV J)	Hc-1	MFIA
Lymphoid Leukosis Viruses	A,B,C,D,E,J	ELISA
Marek's Disease (Serotypes 1,2, 3)	SB-1	AGP
Mycoplasma gallisepticum	A5969	SPA
Mycoplasma synoviae	WVU 1853	SPA
Newcastle Disease	LaSota	MFIA
Reticuloendotheliosis Virus	ATCC 770 T	IFA, AGP
Salmonella pullorum-gallinarum	K Polyvalent	SPA
Salmonella species (Every 4 weeks)		IA

\* Most Recent 1.25% Sampling Date:

No evidence of the listed agents have been detected since the last negative sampling date.

Previous Negative Sampling Dates:

SAMPLE

AGP	Agar Gel Precipitin
CO	Clinical Observation
ELISA	Enzyme-Linked Immunosorbent Assay
HI	Hemagglutination Inhibition
IA	Isolation of Agent
IFA	Indirect Fluorescent Antibody
MNT	Microneutralization
PM	Postmortem
SPA	Serum Plate Agglutination
MFIA	Multiplexed Fluorometric Immunoassay

\* 50 Samples tested per Ph.Eur. 5.2.2

\*\*200 Samples tested per Ph.Eur.5.2.2

taken between the ages of 12-16 weeks and 16-20 weeks.

\_\_\_\_\_  
 Technical Services Supervisor Date  
 Avian Vaccine Services

\_\_\_\_\_  
 Quality Assurance Date  
 Avian Vaccine Services