



AutoClip®

Wound Closure System

The AutoClip Wound Closure System

provides a fast and effective alternative to sutures for closing wounds. It is the perfect system for quickly closing incisions made for ALZET pump implantations, enabling researchers to streamline surgeries and save valuable time. The AutoClip system consists of clip Appliers, Remover, and 9 mm, stainless steel wound clips. The automatic, spring-loaded Applier is easy to use, and the clips provide maximum holding force with minimal animal trauma. When the incision is healed, the wound clips are easily removed using the AutoClip Remover. The 9 mm AutoClips are non-pyrogenic and optimum for use in mice, rats and other laboratory animals. All AutoClip components are supplied non-sterile, but may be autoclaved.



AutoClip Advantages

- ☀ Easy to use
- ☀ Fast and efficient wound closure
- ☀ Automatic, unassisted application
- ☀ Reduces incidence of infection
- ☀ Autoclavable

➔ Order today by contacting our Customer Service at: info@eu.crl.com

Caution: For use in laboratory research animals (not for use in humans)

Item	Description
AutoClip Kits 	A convenient AutoClip kit is available containing the stainless steel Applier, Remover, and 100 wounds clips (Kit 100) or 500 wound clips (Kit 500). Individual components may also be ordered separately. AutoClip Kit 100 (contains 100 clips) Order No.: 0009953 AutoClip Kit 500 (contains 500 clips) Order No.: 0009954
AutoClip Applier 	The AutoClip Applier is an easy to use, stainless steel instrument for fast and efficient application of the 9 mm wound clips. It essentially works like a staple gun, with an automatic, spring-loaded design that enables quick and easy closure of surgical incisions. Order No.: 0009956
AutoClip Remover 	The AutoClip Remover enables quick and easy removal of the wound clips after the incision has healed. Order No.: 0009957
AutoClip Wound Clips 	The 9 mm AutoClips are designed to be used with the AutoClip Applier, and they provide maximum holding force with minimal tissue trauma. They are packaged 20 clips on a cartridge. AutoClips are available in quantities of 100, 500, and 1,000. AutoClips (Qty. 100) Order No.: 0009950 (pkg of 100) AutoClips (Qty. 500) Order No.: 0009951 (pkg of 500) AutoClips (Qty. 1,000) Order No.: 0009952 (pkg of 1,000)