



LAL ACCESSORY PRODUCTS

One of the most common reasons for inconsistent LAL results is contamination and interference associated with uncontrolled accessories. Charles River supplies all the necessary accessory products required to run an LAL test, shipped with Certificates of Quality certifying appropriate endotoxin levels.

LAL Accessory Products

Control Standard Endotoxin (CSE) *E. coli*

	Pack	Code
500 ng per vial*	6/pack	E110
10 ng per vial*	6/pack	E120
Positive product controls (for single test)	25/pack	PC100
Extended CSE Dilution Kit (CXE)*		E140

*Lot-specific Certificate of Analysis included.

LAL Reagent Water

	Pack of 6 Code	Case of 12 Code
30 mL bottle <0.005 EU/mL		W130
50 mL bottle <0.005 EU/mL		W120
100 mL bottle <0.005 EU/mL		W110
500 mL bottle <0.005 EU/mL	W150	

2007

LAL Buffers

	Pack/Case	Code
30 mL 0.1 M Tris Buffer [†]	12/case	BT103
5 mL 0.25 M Tris Buffer [†]	6/pack	BT101
5.5 mL Endotoxin-Specific Buffer [†]	6/pack	BG120

[†]Please contact Technical Assistance before using buffers. Additional buffers are available for specific testing needs.

Endotoxin Indicators**

	Case of 12 Code	Each Code
2,000 EU	EVV2K	
10,000 EU	EVV10K	
100,000 EU	EVV100K	
1 million EU		EVV1M
2.5 million EU		EVV2.5M
10 million EU		EVV10M

**For dry-heat oven validations.



LAL Accessory Products

2007

Endotoxin-Free Test Tubes

	Pack	Code
Depyrogenated flint glass tubes with caps 10 x 75 mm	25/pack	T100
Depyrogenated flint glass tubes wrapped in aluminum foil 10 x 75 mm	50/pack	T200
Depyrogenated borosilicate dilution tubes wrapped in aluminum foil 8 x 75 mm	50/pack	T500
10 x 75 mm	50/pack	T400
13 x 100 mm	32/pack	T300
18 x 150 mm	14/pack	T600

Pipettes and Pipette Tips

	Pack	Code
Depyrogenated glass pipettes*		
1 mL	10/pack	P100
2 mL	10/pack	P200
5 mL	10/pack	P500
10 mL	10/pack	P1000
Eppendorf® pipette tips (individually wrapped)	50/pack	D100*
Eppendorf® pipette tips (individually wrapped)	100/pack	D200

*No certificate of quality available.

96-Well Plates*

	Pack	Code
96-well polystyrene plate TCT (certified to 0.05 EU/mL)	each	M905
96-well polystyrene plate TCT (certified to 0.05 EU/mL)	50/case	M905
96-well polystyrene plate TCT (certified to 0.005 EU/mL)	each	M9005
96-well polystyrene plate TCT (certified to 0.005 EU/mL)	50/case	M9005

*Tissue culture treated.



Charles River offers a full line of accessory products that are certified to be endotoxin-free to support trouble-free testing.