

1.0 DESCRIPTION

Golden brown granular wood product manufactured by first drying cross cut Scottish Pine using a drum blender at approximately 80°C*. The dry material is then sifted to produce a free flowing dust free wood grain.

*Temperature may vary according to the moisture level.

2.0 CHEMICAL PROPERTIES

None of the wood used to produce Datesand Grade 7 has been chemically treated in any way.

3.0 PHYSICAL PROPERTIES

Moisture Content : <12.0 % Post Production
Loose Bulk Density : 110 - 210 grams/litre

Particle Size Using Sieve Shaker (WTI803QC)

<u>Mesh</u>	<u>Microns</u>	<u>Cumulative % Retained</u>	
		<u>Minimum</u>	<u>Maximum</u>
6	2800	2.0	10.0
16	1000	75.0	98.0

4.0 PACKAGING

Packed in 110 litre approx. 10kg autoclavable polypropylene sacks marked with Product Type/Grade, Batch number and Volume. Also available in 180kg tote sacks.