

LABORATORY GRADE ASPEN and SHREDDED ASPEN SHAVINGS

Contact Bedding Processed from
Virgin Aspen Hardwood



- Compatible with all animals, including transgenic breeds
- Heat-treated to reduce bacteria and enhance absorbency
- Proven track record of superior performance
- Shredded variety is particularly suited for nesting and animal enrichment
- Aspirated to remove dust
- Regular changing prevents the development of ammonia
- Creates a comfortable habitat with a natural appearance
- Proven performance for micro-isolator cage systems



Northeastern Products Corp.
115 Sweet Road • P.O. Box 98
Warrensburg, NY 12885

T 518.623.3161
800.873.8233
F 518.623.3803
E info@labbedding.com
W labbedding.com

Premium Lab-Bedding Since 1955!

Natural
Bedding
Solutions . . .
for
Today's
Research!

LABORATORY GRADE ASPEN and SHREDDED ASPEN SHAVINGS

Contact Bedding Processed from
Virgin Aspen Hardwood



Lab Grade



Shredded

Raw Material — Virgin, Aspen
Hardwood

Standard Particle Size

- Less than 5% retained on a
1 inch x 1 inch screen
- Less than 3% passes a #20
ASTM Sieve

Moisture Content — $9 \pm 3\%$

Package Size

- 8.0 Cubic Feet loose,
compressed to 2.6 cubic feet

Package Type

- Two-ply autoclavable paper bag
- All bags are production date
coded for inventory control

Laboratory Analysis

- Performed annually; results are
available upon request
- Consistently meets or exceeds
industry standard contaminant
specifications for contact bedding

Call 800.873.8233 for details
and laboratory analysis

**Laboratory Beddings Available
from Northeastern Products Corp.:**

- Alpha Chip[®]
- Aspen Chip
- Beta Chip[®]
- Aspen Shavings
- Pine Shavings
- Corn Cob
- SoftZorb[®]
- Certified Bedding
- Irradiated Bedding

Premium Lab-Bedding Since 1955!

