

Miscellaneous Antiserum/Antibody

Rabbit Antiserum

Rabbit p27 antiserum was prepared in rabbits using purified virion protein (m.w.=27,000) extracted from avian myeloblastosis virus. Antiserum reacts with all subgroups of Leukosis virus.

Rabbit p27 Antiserum	Material #
Liquid	10100766
Lyophilized	10100763

These reagents are derived from our Rabbit p27 antiserum. Normal useful working dilution of Anti p27 IgG and Anti p27 conjugate are approximately 1:1000

Rabbit p27 Antiserum	Material #
Rabbit Anti-p27 IgG (Lyophilized)	10100771
Rabbit Anti-p27 IgG Conjugated to HRP (Lyophilized)**	10100770
p27 Positive TC Control (Lyophilized)	10100764
p27 Negative TC Control (Lyophilized)	10100765

** Restoring Buffer (30% glycerin) supplied at no additional charge when purchasing conjugate. Material # 10100705 if sold separately.